

FORE-SIGHT[®] Roll Stand

Roll Stand with Basket and Padded Grip Handle

CAS Catalog Number: 01-06-0238

Specifications:

- Load capacity: 30 lbs (13.6 kg)
- Overall height to mounting plate: 43" (109 cm)
- 21" (53.3 cm) diameter aluminum base
- 4" (10 cm) dual wheel casters (2 locking, 3 non-locking) for easy transition over uneven surfaces
- Built-in 10 lbs (4.5 kg) counterweight (beneath the base) for tip safety
- 5" (12.7 cm) Slide – In style mounting plate
- Adjustable tilt angle for optimum viewing
- Utility basket with hooks for cable management
12" w x 6" h x 8" d (30.5 x 20.3 x 15.2 cm)
- Padded grip handle
- Requires FORE-SIGHT adapter plate:
P/N 01-06-0237





01-06-0237
Adapter Plate



01-06-0223
FS 6" Utility Basket
For attachment to 16" & 8" Arm wall mount
and 16" Arm pole mount



01-06-0232
FS 3" Utility Basket
For attachment to 16" & 8" Arm wall mount
and 16" Arm pole mount



01-06-0236
FS 16" Arm Pole Mount
Adapter Plate to be purchased separately



01-06-0235
FS 8" Arm Wall Mount
(Giraffe Bed)
Adapter Plate to be purchased separately



01-06-0234
FS 16" Arm Wall Mount
Adapter Plate to be purchased separately



01-06-0231
Giraffe Dovetail Attachment



*** NEW - Special order**
FS 8" Arm Pole Mount
Adapter Plate to be purchased separately



28-02-0555
8" M-Series w/ 5" P-Plate
for 2" Post